MODEL 253 AND 253-L (WATERMARK 200) SOIL MOISTURE SENSOR

REVISION: 6/96

COPYRIGHT (c) 1993-1996 CAMPBELL SCIENTIFIC, INC.

WARRANTY AND ASSISTANCE

The **MODEL 253 AND 253-L (WATERMARK 200) SOIL MOISTURE SENSOR** is warranted by CAMPBELL SCIENTIFIC, INC. to be free from defects in materials and workmanship under normal use and service for twelve (12) months from date of shipment unless specified otherwise. Batteries have no warranty. CAMPBELL SCIENTIFIC, INC.'s obligation under this warranty is limited to repairing or replacing (at CAMPBELL SCIENTIFIC, INC.'s option) defective products. The customer shall assume all costs of removing, reinstalling, and shipping defective products to CAMPBELL SCIENTIFIC, INC. CAMPBELL SCIENTIFIC, INC. will return such products by surface carrier prepaid. This warranty shall not apply to any CAMPBELL SCIENTIFIC, INC. products which have been subjected to modification, misuse, neglect, accidents of nature, or shipping damage. This warranty is in lieu of all other warranties, expressed or implied, including warranties of merchantability or fitness for a particular purpose. CAMPBELL SCIENTIFIC, INC. is not liable for special, indirect, incidental, or consequential damages.

Products may not be returned without prior authorization. To obtain a Returned Materials Authorization (RMA), contact CAMPBELL SCIENTIFIC, INC., phone (435) 753-2342. After an applications engineer determines the nature of the problem, an RMA number will be issued. Please write this number clearly on the outside of the shipping container. CAMPBELL SCIENTIFIC's shipping address is:

CAMPBELL SCIENTIFIC, INC. RMA#_____ 815 West 1800 North Logan, Utah 84321-1784

CAMPBELL SCIENTIFIC, INC. does not accept collect calls.

Non-warranty products returned for repair should be accompanied by a purchase order to cover the repair.



CAMPBELL SCIENTIFIC, INC.

815 W. 1800 N. Logan, UT 84321-1784 USA

Phone (435) 753-2342 FAX (435) 750-9540 www.campbellsci.com Campbell Scientific Canada Corp. 11564 -149th Street Edmonton, Alberta T5M 1W7 CANADA Phone (780) 454-2505 FAX (780) 454-2655 Campbell Scientific Ltd. Campbell Park 80 Hathern Road Shepshed, Loughborough LE12 9GX, U.K. Phone +44 (0) 1509 601191

MODEL 253 AND 253-L (WATERMARK 200) SOIL MOISTURE SENSOR

1. GENERAL

The Watermark 200 (CSI sensor Models 253, 253-L, 257, and 257-L) provides a convenient method of estimating water potential between 0 and 2 bars (wetter soils) with a Campbell Scientific CR10, 21X, or CR7 datalogger. CSI Models 253 and 253-L are for connection to the AM32 or AM416 Analog Multiplexers. Models 257 and 257-L connect directly to a datalogger.

The Watermark block estimates water potential. For applications requiring high accuracy, call a Campbell Scientific applications engineer for information on precision soil moisture measurement systems.

The Watermark consists of two concentric electrodes embedded in a reference matrix material. The matrix material is surrounded by a synthetic membrane for protection against deterioration. An internal gypsum tablet buffers against the salinity levels found in irrigated soils.

If cultivation practices allow, the sensor can be left in the soil all year, eliminating the need to remove the sensor during the winter months.



FIGURE 1.1 253 Soil Moisture Sensor

2. INSTALLATION AND REMOVAL

Placement of the Watermark is important. To acquire representative measurements, avoid high spots, slope changes, or depressions where water puddles. Typically, the sensor must be located in the root system of the crop.

- Soak the sensors overnight in irrigation water. Always install a wet sensor. If time permits, allow the sensor to dry for 1 to 2 days after soaking, and repeat the soak/dry cycle twice to improve sensor response.
- 2. Make a sensor access hole to the depth required with a 7/8" rod. Fill the hole with water and push the sensor to the bottom of the hole. Very coarse or gravelly soils may require an oversized hole (1 to 1-1/4") to prevent abrasion damage to the sensor membrane. In this case, you will need to "grout in" the sensor with a slurry made from the sample soil to get a snug fit in the soil.

Snug fit in the soil is most important. Lack of a snug fit is the premier problem in sensor effectiveness. In gravelly soils, and with deeper sensors, sometimes it is hard to get the sensor in without damaging the membrane. The ideal method of making the access hole is to have a "stepped" tool that makes an oversized hole for the upper portion and an exact size hole for the lower portion. In either case, the hole needs to be carefully backfilled and tamped down to prevent air pockets which could allow water to channel down to the sensor.

A length of 1/2" class 315 PVC pipe fits snugly over the sensor collar and can be used to push in the sensor.

You can leave the PVC in place with the wires threaded through the pipe and the open end taped shut (duct tape is adequate). This practice also makes it easy to remove sensors used in annual crops. When doing this, solvent weld the PVC pipe to the sensor collar. Use PVC/ABS cement on the stainless steel sensors with the green top. Use clear PVC cement only on the PVC sensors with the gray top.

- 3. When removing sensors prior to harvest in annual crops, do so just after the last irrigation when the soil is moist. Do not pull the sensor out by the wires. Careful removal prevents sensor and membrane damage.
- 4. When sensors are removed for winter storage, clean, dry, and place them in a plastic bag.

NOTE: The black outer jacket of the cable is Santoprene[®] rubber. This compound was chosen for its resistance to temperature extremes, moisture, and UV degradation. However, this jacket will support combustion in air. It is rated as slow burning when tested according to U.L. 94 H.B. and will pass FMVSS302. Local fire codes may preclude its use inside buildings.

3. WIRING

The model 253 sensor is supplied with two green leads from Watermark. The leads from the Watermark electrode are connected directly to the H and L inputs on the AM32 or AM416. The lead coming from the center of the sensor is connected to H and the lead from the outer portion of the sensor to L. The wires can be differentiated by the grooved strip in one of the leads of the green wires. On the 253-L, Campbell Scientific splices a two conductor shielded cable to the two conductor green cable supplied from Watermark. The black conductor is connected to H, the white conductor to L, and the shield wire to G or $\frac{1}{2}$. A 1k Ω resistor at the datalogger is used to complete the half bridge measurement.

4. MEASUREMENT

Instruction 5, AC Half Bridge, is used to excite and measure the model 253. Recommended excitation voltages and input ranges are listed in Table 1.

TABLE 1. Excitation and Voltage Range

<u>DATALOGGER</u>	<u>mV EX</u>	RANGE <u>CODE</u>	<u>FSR</u>
21X	500	14	± 500 mV
CR10	250	14	± 250 mV

NOTE: Do not use a slow integration time as sensor polarization errors will occur.

4.1 CALCULATE SENSOR RESISTANCE -INSTRUCTION 59

Instruction 59, Bridge Transform, is used to output sensor resistance (R_s). The instruction takes the AC Half Bridge output (Vs/Vx) and computes the sensor resistance as follows.

$$R_{s} = R1(X/(1-X))$$

where, X = Vs/Vx (Output from Instruction 5)

A multiplier of 1 should be used to output sensor resistance (R_s) in terms of $k\Omega$.

4.2 CALCULATE SOIL WATER POTENTIAL

The datalogger can calculate soil water potential (bars) from the sensor resistance (R_s) and soil temperature (T_s). See Table 2.

The need for a precise soil temperature measurement should not be over emphasized. Soil temperatures vary widely where placement is shallow and solar radiation impinges on the soil surface. A soil temperature measurement may be needed in such situations, particularly in research applications. Many applications, however, require deep placement (5 to 10 inches) in soils shaded by a crop canopy. A common practice is to assume the air temperature at sunrise will be close to what the soil temperature will be for the day.

4.2.1 Linear Relationship

For applications in the range of 0 to 2 bars, the water potential and temperature responses of the Watermark can be assumed to be linear (measurements beyond 1.25 bars have not been verified, but work in practice).

The following equation normalizes the resistance measurement to 21°C.

$$R_{21} = \frac{R_s}{1 - (0.018 * dT)}$$
 [1]

where

 R_{21} = resistance at 21°C

R_s = the measured resistance

T_s = soil temperature

Water potential is then calculated from R_{21} with the relationship.

SWP =
$$0.07407 * R_{21} - 0.03704$$
 [2]

SWP = Soil Water Potential (bars)

4.2.2 Non-Linear Relationship02For more precise work, calibration and temperature
compensation in the range of 0.1 to 1.00 bar has
been refined by Thompson and Armstrong (1987),
as defined in the non-linear equation,03SWP = $\frac{R_s}{0.01306[1.062(34.21 - T_s + 0.01060T_s^2) - R_s]} *.01$ [3]Table 2. Comparison of Estimated Soil
Water Potential and Rs at 21°C04

Bars (Non- Linear <u>Equation)</u>	Bars (Linear <u>Equations)</u>	<u>(Rs)kOhms</u>
.09 .14 .20 .27 .35 .45 .56 .69 .85 1.05	.037 .11 .18 .26 .33 .41 .48 .56 .63 .70 .78 .85 .92 .99 1.07 1.15 1.22 1.29 1.44 1.59 1.74 1.88 1.99	1.00 2.00 3.00 4.00 5.00 6.00 7.00 8.00 9.00 10.00 11.00 12.00 13.00 14.00 15.00 16.00 17.00 18.00 20.00 22.00 24.00 26.00 27.50

02: 01: 02:	P87 0 32	Beginning of Loop Delay Loop Count
03: 01: 02: 3] 03: 04:	P22 1 1 5000	Excitation with delay (clock) Excitation Channel #1 Excite for 0.01 seconds 0 second delay after excitation Excitation = 5000 mV
04:	P5	AC Half Bridge (Measure AC Conductivity)
01: 02: 03: 04:	1 14 1 2	Rep 500 mV Fast Range In Channel Excite All Reps w/Excite Channel 2
05: 06:	500 1	mV excitation Location (Indexed Location to Store) [:kOhms#1]
07: 08:	1 0	Multiplier Offset
05:	P95	End
06: 01: 02:	P20 0 1	Set Port (Reset AM32) Set Low Port Number
07:	P59	BR Transform Rf[X/(1-X)] (Compute Resistances)
01: 02: 03:	32 1 1	Reps Location [:kOhms#1] Multiplier (Rf/1000)




FIGURE 5.1

5. PROGRAMMING (MEASURING BLOCK RESISTANCE)

The following examples demonstrate the connections and programming used to measure the resistances (kohms) of 32 soil moisture blocks.

5.1 AM32 AND 21X

See Figure 5.1 for wiring diagram.

01:	P20	Set Port (Enable AM32)
01:	1	Set High `
02:	1	Port One

5.2 AM32 AND CR10

See Figure 5.2 for wiring diagram. This program can also be used with 21Xs with OSX PROMs although the clock pulse delay is 0.1 seconds (Figure 5.3).

02:P87Beginning of Loop01:0Delay02:32Loop Count	
03: P86 Do (Clock Pulse, 10 ms) 01: 76 Pulse Port 6	
04: P5 AC Half Bridge (Measure AC Conductivity)	
01: 1 Rep	
02: 14 250 mV Fast Range	
03: 1 In Channel	
04: 2 Excite All Reps w/ Excitation Channel 2	
05: 250 mV Excitation	
06: 1 Location (Indexed Location	n
to Store) [:kOhms#1]	
07: 1 Multiplier	
08: 0 Offset	
05: P95 End	
06: P86 Do (Reset AM32)	
01: 55 Set Port 5 low	
07: P59 BR Transform Rf[X/(1-X)] (Compute Resistances)	
01: 32 Reps	
02: 1 Location [:kOhms#1]	
03: 1 Multiplier (Rf/1000)	






FIGURE 5.2



5.3 AM416 AND 21X

See Figure 5.4 for wiring diagram.

qq
elay (clock) nel #1 econds after excitation 0 mV
(Measure
nge v/
nel 2 d Location s#1]
AM416)
f[X/1(1-X)] ances s#1] 00)






5.4 AM416 AND CR10

See Figure 5.5 for wiring diagram. This program can also be used with 21X with OSX PROMs although the clock pulse delay is 0.1 seconds (Figure 5.6).

01: 01:	P86 41	Do Set Port (Enable AM416) Set High Port 1
02: 01: 02:	P87 0 16	Beginning of Loop Delay Loop Count
03: 01:	P86 72	Do (Clock Pulse, 10 ms) Pulse Port 2
04: 01:	P90 2	Step Loop Index Step
05:	P5	AC Half Bridge (Measure AC Conductivity)
01:	2	Rep
02:	14	250 mV Fast Range
03:	1	In Channel
04:	2	Excite All Reps w/ Excitation Channel 2
05:	250	mV Excitation
06:	1	Location (Indexed Location to Store) [:kOhms#1]
07:	1	Multiplier
08:	0	Offset
06:	P95	End
07:	P86	Do (Reset AM416)
01:	51	Set Low Port 1
08: 01: 02:	P59 32 1	BR Transform Rf[X/(1-X)] Reps Location [:kOhms#1]
03:	1	Multiplier (Rf/1000)








Figure 5.6

6. PROGRAMMING (CALCULATING SOIL WATER POTENTIAL)

6.1 LINEAR RESISTANCE AND TEMPERATURE RELATIONSHIP (0 TO 2 BARS)

Calculate Temperature Correction Factor. See Equation [1] in Section 4.2.1...

...Calculate dT = T - 21

04:	P34	Z=X+F
01:	34	X Loc TmpDegC
02:	-21	F
03:	36	Z Loc [:CorrFactr]
Calcu	ulate (0.018	* dT)
05:	P37	Z=X*F
05: 01:	P37 36	Z=X*F X Loc CorrFactr
		= /
01:	36	X Loc CorrFactr

253 AND 253-L SOIL MOISTURE SENSOR

	06: 01:	P34 36	Z=X+F X Loc CorrFactr	07:	P30	Z=F
	01:	-1	F	07:	34.21	F
	03:	36	Z Loc [:CorrFactr]	02:	35	Z Loc [:Constant]
	07:	P37	Z=X*F		04 04 T-	
	01:	36	X Loc CorrFactr	5VVP =	= 34.21-Ts	
	02:	-1	F	08:	P35	Z=X-Y
	03:	36	Z Loc [:CorrFactr]	01:	35	X Loc Constant
	Apply T	emperature	Correction and	02:	34	Y Loc Tsoil C
			to Ohm Measurements. See	03:	35	Z Loc [:Constant]
	Equatio	n [2] in Sect	ion 4.2.1			T-) · /T- 40*0 04000)]
	08:	P87	Beginning of Loop	SVVPC	alc.=[(34.21-	Ts)+(Ts^2*0.01060)]
	01:	0	Delay	09:	P33	Z=X+Y
	02:	32	Loop Count	01:	41	X Loc Bars#1
	Temp	erature Cori	rect Ohms:	02:	3	Y Loc Constant
	09:	P38	Z=X/Y	03:	41	Z Loc [:Bars#1]
	01:	1	X Loc kOhms#1	S\//D		
	02:	36	Y Loc CorrFactr	3007 =	= 1.062[SWP	calc.j
	03:	41	Z Loc [:Bar#1]	10:	P37	Z=X*F
	Apply	Calibration	Slope and Offset	01:	41	X Loc Bars#1
	10:	P37	Z=X*F	02:	1.062	F
	01:	41	X Loc Bar#1	03:	41	Z Loc [:Bars#1]
	02:	.07407		11:	Р	35 Z=X-Y
	03:	41	Z Loc [:Bar#1]		1	55 Z=X 1
	11:	P34	Z=X+F	SWP =	= SWPcalcF	Rs
	01:	41	X Loc Bar#1			
	02:	03704		01:	41	X Loc Bars#1
	03:	41	Z Loc [:Bar#1]	02: 03:	1 41	Y Loc kOhms#1 Z Loc [:Bars#1]
	12:	P95	End	00.	71	2 LOC [.Da13#1]
				SWP =	= 0.01306*SV	VPcalc.
6.2						
	BAR)	RAIURER	ELATIONSHIP (0.1 TO 1	12:	P37	Z=X*F
				01: 02:	41 0.0130	X Loc Bars#1 F
			ctions convert R _S to Soil	02:	41	Z Loc [:Bars#1]
	Water F Section		pars. See Equation [3] in			[]
				SWP =	= Rs/SWPca	c.
	08:	P87	Beginning of Loop	40.	DOO	7
	01: 02:	0 32	Delay Loop Count	13: 01:	P38 1	Z=X/Y X Loc kOhms#1
	02.	32	Loop Count	01.	41	Y Loc Bars#1
	05:	P36	Z=X*Y SWP = Tsoil^2	03:	41	Z Loc [:Bars#1]
	01:	34	X Loc Tsoil C			
	02:	34	Y Loc Tsoil C	SWPb	ars = SWP(k	:Pa) * 0.01
	03:	41	Z Loc [:Bars#1]	14:	P37	Z=X*F
	SWP =	Tsoil^2 * 0.0	0106	01:	41	X Loc Bars#1
				02:	0.01	F
	06: 01:	P37 41	Z=X*F X Loc Bars#1	03:	41	Z Loc [:Bars#1]
	01.	0.0106	F		D 0-	
	03:	41	Z Loc [:Bars#1]	15:	P95	End Calculation Loop
	•					

253 AND 253-L SOIL MOISTURE SENSOR

7. PROGRAMMING (COMPREHENSIVE)

Follow these steps to create a complete program:

- <u>Step 1</u>. Allocate at least 75 input locations in EDLOG.
- <u>Step 2</u>. Set the execution interval according to need:
- * 1 Table 1 Programs 01: 3600 Sec. Execution Interval
- <u>Step 3</u>. Make a temperature measurement to correct for temperature effects. Select from options 1 or 2.
 - <u>Option 1</u>. If a 107B Probe is part of your system, measure soil temperature:

01:	P11	Temp 107 Probe
01:	1	Rep
02:	3	IN Chan
03:	3	Excite all reps w/EXchan 3
04:	34	Loc [:TempDegC]
05:	1	Mult
06:	0	Offset

Option 2. If a 107B Probe is not available but a 107 Probe is part of your system, measure air temperature in the early morning (6:00 A.M.) and assume that will be the soil temperature for the day:

02: 01: 02: 03:	P92 360 1440 30	If time is minutes (seconds) into a minute or second interval Then Do
03: 01: 02: 03: 04: 05: 06:	P11 1 3 34 1 0	Temp 107 Probe Rep IN Chan Excite all reps w/EXchan 3 Loc [:TempDegC] Mult Offset
04:	P95	End

<u>Step 4</u>. Make the resistance measurements. It may be appropriate to make measurements only once or twice a day. This example makes measurements twice a day at 6:00 AM and 6:00 PM:

05:	P92	If time is
01:	360	minutes (seconds) into a
02:	720	minute or second interval
03:	30	Then Do

Insert one of the examples from Section 5 here.

<u>Step 5</u>. Calculate soil water potential using resistance and temperature:

Insert one of the examples from Section 6 here.

Step 6. Output data to final storage after each measurement and calculation:

06:	P86	Do
01:	10	Set high Flag 0 (output)
07:	P77	Real Time
01:	0220	Day,Hour-Minute
08:	P70	Sample kOhm Resistances
01:	32	Reps
02:	1	Loc
09:	P70	Sample Deg C Temperature
01:	1	Reps
02:	34	Loc
10:	P70	Sample Bar Potential
01:	32	Reps
02:	41	Loc
11:	P95	End 6 am and 6 pm loop

8. INTERPRETING RESULTS

As a general guide, Watermark 200 measurements indicate soil moisture as follows:

- 0 to 10 centibars = Saturated soil.
- 10 to 20 centibars = Soil is adequately wet (except coarse sands, which are beginning to lose water).
- 30 to 60 centibars = Usual range for irrigation (except heavy clay).
- 60 to 100 centibars = Usual range for irrigation for heavy clay soils.
- 100 to 200 centibars = Soil is becoming dangerously dry for maximum production.

9. TROUBLESHOOTING

To test the sensor, submerge it in water. Measurements should be from -.03 to .03 bars. Let the sensor dry for 30 to 48 hours. You should see the reading increase from 0 to 150+. Put the sensor back in the water. The reading should run right back down to zero in 1 to 2 minutes. If the sensor passes these tests, consider the following.

- Sensor may not have a snug fit in the soil. This usually happens when an oversized access hole has been used and the backfilling of the area around the sensor is not complete.
- Sensor is not in an active portion of the root system, or the irrigation is not reaching the sensor area. This can happen if the sensor is sitting on top of a rock or below a hard pan which may impede water movement. Re-installing the sensor usually solves this problem.

3. When the soil dries out to the point where you are seeing readings higher than 80 centibars, the contact between soil and sensor can be lost because the soil may start to shrink away from the sensor. An irrigation which only results in a partial rewetting of the soil will not fully rewet the sensor, which can result in continued high readings from the Watermark. Full rewetting of the soil and sensor usually restores soil/sensor contact. This is most often seen in the heavier soils and during peak crop water demand when irrigation may not be fully adequate. The plotting of readings on a chart is most useful in getting a good picture of this sort of behavior.

Reference

Thompson, S.J. and C.F. Armstrong, Calibration of the Watermark Model 200 Soil Moisture Sensor, Applied Engineering in Agriculture, Vol. 3, No. 2, pp. 186-189, 1987.

Parts of this manual were contributed by Irrometer Company, Inc., manufacturer of the Watermark 200.